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Creating the future together.

CONTENTS

1	SAFETY PRECAUTIONS	1
2	GENERAL INFORMATION	2
2.1	Description	2
2.2	Safety precautions	2
2.3	Technical data	2
2.4	Accessories supplied with centrifuge	3
2.5	Warranty	3
3	INSTALLATION	3
3.1	Unpacking the centrifuge	3
3.2	Required space	3
3.3	Installation	3
4	INSTALLATION OF ROTORS AND ROTOR MAINTENANCE	4
4.1	Rotors and accessories	4
4.2	Rotor maintenance	5
4.3	Removing and installing the angle rotor	5
4.4	Loading the rotor	5
4.5	Overloading rotors	5
5	OPERATION	6
5.1	Attaching rotor lid	6
5.2	Closing the lid	6
5.3	Lid release	6
5.4	Lid lock	6
5.5	Speed selection (see Figure 2)	6
5.6	Selection of operating time, momentary operation, Start/Stop	7
6	SERVICE AND MAINTENANCE	7
6.1	Centrifuge service	7
6.2	Cleaning the centrifuge	7
6.3	Cleaning the rotor	7
6.4	Disinfection of the rotor	7
6.5	Replacing fuses	7
7	TROUBLESHOOTING	8
8	TECHNICAL SUPPORT	9
9	DETERMINATION OF G-VALUES	9

1 SAFETY PRECAUTIONS

- NEVER** use the centrifuge in any manner not specified in these instructions.
- NEVER** operate the centrifuge without a rotor properly attached to the shaft.
- NEVER** tighten rotor nut by hand only.

- NEVER** fill tubes while they are in the rotor. Liquid spillage may harm unit.
- NEVER** put hands in the rotor area unless the rotor is completely stopped.

- NEVER** move the centrifuge while the rotor is spinning.

- NEVER** use solvents or flammables near this or other electrical equipment.
- NEVER** centrifuge flammable, explosive or corrosive materials.

- NEVER** centrifuge hazardous materials outside of a hood or proper containment facility.

- ALWAYS** load the rotor symmetrically. Each tube should be counter balanced by another tube of the same type and weight.

- ALWAYS** locate the centrifuge within easy access to an electrical outlet.

- ALWAYS** use only microcentrifuge tubes made from plastic and designed to withstand centrifugal forces of at least 21,200 x g.

- ALWAYS** use a wrench to tighten the rotor nut.

2 GENERAL INFORMATION

This manual provides important safety information for the **PerfectSpin 24 Plus** laboratory microcentrifuge (Cat. No. 91-PSPIN-24P/91-PSPIN24P-US). It should be kept near the centrifuge for quick and easy reference.

2.1 Description

The **PerfectSpin 24 Plus** is a small benchtop centrifuge designed for separation of various research samples. The motor is brushless and requires no routine maintenance. The **PerfectSpin 24 Plus** is supplied with a 24 x 1.5 / 2.0 ml rotor for micro samples. Adapters are available for tubes smaller than 1.5 ml. The **PerfectSpin 24 Plus** reaches speeds of up to 15,000 rpm / 21,200 x g.

2.2 Safety precautions

Note: All users of the centrifuge must read the Safety Precautions section of this manual before attempting to operate the unit!

If this equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Do not operate the centrifuge if any of the following conditions exist:

- The centrifuge has not been installed properly.
- The centrifuge is partially dismantled.
- Service has been attempted by unauthorized or unqualified personnel.
- The rotor has not been installed securely on the motor shaft.
- Rotors and accessories not belonging to the standard range are being used without permission being obtained from the manufacturer to use such rotors and / or accessories in the centrifuge. Exception: Microcentrifuge tubes made of plastic, normally available in the laboratory.
- The centrifuge is located in an explosive atmosphere.
- Materials to be centrifuged are combustible and/or explosive.
- Materials to be centrifuged are chemically reactive.
- The rotor load is not properly balanced.
- The rotor nut was not tightened with a wrench.

2.3 Technical data

Dimensions:

Width	9.45 inches
Depth	13.8 inches
Height	7.48 inches
Maximum speed	15,000 rpm
Maximum RCF	21,200 x g
Maximum volume	24 x 1.5 / 2.0 ml
Admissible sample density	1.2 kg / dm ³
Electrical/fuse rating	120 V~, 50-60 Hz, 1.9 A, 5 AT 230 V~, 50-60 Hz, 1.1 A, 2.5 AT
Operating temp./humidity	0 °C to 40 °C – ≤ 80 % RH

2.4 Accessories supplied with centrifuge

Each unit is supplied with 1 instruction manual, 1 power cord, a standard rotor and a rotor removal tool.

2.5 Warranty

PEQLAB guarantees that the **PerfectSpin 24 Plus** laboratory microcentrifuge you have received has been thoroughly tested and meets its published specification.

However, immediately upon arrival, please check carefully that the shipment is complete and has not been damaged in transit. For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SUPPORT'). Please retain all packaging materials until the delivery has been completely checked since this will speed up the return of goods if required and reduce environmental impact. Any form of returns, replacements or credit notes must be agreed in advance by PEQLAB.

This centrifuge has been subjected to thorough testing and quality control. In the unlikely event of a manufacturing fault, our two year warranty (from the date of delivery) covers the centrifuge and the rotor. This warranty becomes invalid in the case of incorrect operation, use of non-standard spare parts or accessories and unauthorized modification of the rotor or centrifuge.

PEQLAB reserves the right to make technical modifications.

3 INSTALLATION

3.1 Unpacking the centrifuge

Before unpacking the centrifuge, inspect the outside of the carton for any shipping damage.

The centrifuge is delivered in a carton with protective foam cushioning. Remove the centrifuge from the carton. Retain the carton and cushioning until it has been established that the centrifuge is working properly.

For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SERVICE').

The accessories supplied with the centrifuge should be kept with the instruction manual near the centrifuge's place of installation.

3.2 Required space

The centrifuge should be installed on a rigid, even surface such as a stable laboratory bench, countertop, etc. To guarantee sufficient ventilation, ensure that the centrifuge has at least 15 cm (6 inches) of free space on all sides, including the rear.

The centrifuge should not be located in areas subject to excessive heat such as in direct sunlight or near radiators or the exhaust of a compressor, as a buildup of heat may occur within the chamber.

3.3 Installation

Before operating the centrifuge, check that the power source (electrical outlet on the wall) corresponds to that on the manufacturer's rating label, then connect the power cord to the centrifuge and the power source.

4 INSTALLATION OF ROTORS AND ROTOR MAINTENANCE

4.1 Rotors and accessories

The following accessories are included or available for the **PerfectSpin 24 Plus**:

Angle rotor for 24 x 1.5 ml tubes (Included)

Tube measurement	1.5 ml / 2.0 ml (10 x 40 mm)
Max. Speed	15,000 rpm
Centrifuging radius	8.4 cm
Max. RCF (g-value)	21,200 x g

Rotor Removal Tool (Included)

Order no.	91-PSPIN-RW
-----------	-------------

Rotor Lid (Included)

Order no.	91-PSPIN-RL
-----------	-------------

Rotor Securing Screw (Included)

Order no.	91-PSPIN-RS
-----------	-------------

Strip Spin Adapter

Order no.	91-PSPIN-AS
Tubes accepted	0.2 ml tubes or 8 x 0.2 ml strips
Max Speed	15,000 rpm
Range Centrifuge radius	4.32 cm to 5.33 cm
Range RCF (g value)	10,874 x g to 13,405 x g

Adapter for 0.5 ml tubes (pk./6)

Order no.	91-PSPIN-A05
Tube measurement	8 x 30 mm
Max. Speed	15,000 rpm
Centrifuging radius	7.53 cm
RCF (g-value)	18,942 x g

Adapter for 0.4 ml tubes (pk./6)

Order no.	91-PSPIN-A04
Tube measurement	6 x 47 mm
Max. Speed	15,000 rpm
Centrifuging radius	8.4 cm
RCF (g-value)	21,200 x g

Adapter for 0.2 ml tubes (pk./6)

Order no.	91-PSPIN-A02
Tube measurement	6 x 21 mm
Max. Speed	15,000 rpm
Centrifuging radius	7.03 cm
RCF (g-value)	17,684 x g

4.2 Rotor maintenance

The rotor should be cleaned thoroughly after each use. **Thorough cleaning must be performed when spinning samples containing phenol or phenol chloroform.** Periodically inspect the rotor for dents, dings, scratches, discoloration and cracks. If any damage to the rotor is found, discontinue use of the rotor immediately and replace.

4.3 Removing and installing the angle rotor

The **PerfectSpin 24 Plus** comes complete with a standard 24 place rotor installed. To remove the rotor for cleaning, remove the rotor securing screw from the motor shaft by turning the screw counterclockwise, using the rotor wrench. Lift the rotor directly upward in a straight vertical motion.

To replace rotor, first make sure the motor shaft and rotor mounting hole are clean. Place the rotor on the motor shaft. Reinstall the rotor securing screw on the motor shaft by turning it clockwise. Hold the rotor with one hand and **tighten the rotor securing screw, using the rotor wrench.**

4.4 Loading the rotor

Tubes to be loaded should be filled equally by eye. The difference in the weight between the tubes should not exceed 0.1 gram. Tubes should always be loaded so that there is equal spacing between all tubes. One or two additional loaded tubes may need to be added to achieve this. Refer to **Figure 1** to see one typical balancing scheme.

Figure 1: Loading the rotor to achieve balance

4.5 Overloading rotors

The maximum load of the rotor and the maximum speed has been established by the manufacturer. Do not attempt to exceed these values. The maximum speed of the rotor has been established for liquids having a homogeneous density of 1.2 g/ml or less. In order to centrifuge liquids with a higher density it is necessary to reduce the speed. **Failure to reduce the speed may result in damage to the rotor and centrifuge.** The revised maximum speed can be calculated with the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{1.2}{\text{higher_density_value}}} \times \text{max speed } (n_{\text{max}})$$

Example:

Where the density of the liquid is 1.7, the new maximum speed would be calculated as follows:

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.7}} \times 15,000 \text{ rpm} = 12,602 \text{ rpm}$$

If in doubt concerning maximum speeds, please contact the manufacturer for assistance.

5 OPERATION

ATTENTION: Never attempt to operate the centrifuge with rotors or adapters that show signs of corrosion or mechanical damage. Never centrifuge strongly corrosive materials that may damage the rotors, accessories, or bowl of the unit.

5.1 Attaching rotor lid

After the rotor has been properly secured and loaded, attach the rotor lid to the rotor. Always use the rotor lid for safety and to allow the rotor to reach proper speed. Make sure that the rotor lid snaps securely into place, by pressing down on center catch.

5.2 Closing the lid

Close the centrifuge lid. The **PerfectSpin 24 Plus** has a lid lock that activates only when a run is started.

5.3 Lid release

The lid will remain locked during a centrifuge run. Once the run has been completed and the rotor has come to a stop, a beep will indicate the end of a run, and the lid will unlock automatically.

WARNING: Do not attempt to open the lid of any centrifuge until the rotor has come to a complete stop.

In the event of a power failure or malfunction, it may be necessary to open the lid manually.

1. Disconnect the power cord from the wall socket.
2. Remove the plastic plug, located on the left side of the unit.
3. Pull the wire (attached to the plug) to open the lid lock manually.

5.4 Lid lock

The centrifuge can be started only with the lid securely closed. When a run is started, the lid lock automatically activates. Do not attempt to open the lid during a centrifuge run. At the end of the run, the lid will automatically unlock.

Never attempt to override the lid lock mechanism. Doing so is dangerous and could damage the centrifuge.

5.5 Speed selection (see Figure 2)

The speed (rpm or g-force) can be selected from 500 to 15,000 rpm in 100 rpm increments or from 100 to 21,200 x g with the control knob. The speed is selected by pressing the RPM or RCF button. The speed signal will begin to blink. Then, turn the control knob to increase or decrease the value.

Figure 2: PerfectSpin 24 Plus control panel layout

5.6 Selection of operating time, momentary operation, Start/Stop

Operating time can be selected from 0.5 min to 99 min by pressing the TIME button and adjusting with the control knob. The time can be set in 0.5 minute increments from 0 to 10 minutes and in 1 minute increments from 10 to 99 minutes. After 99 minutes, the display shows "--" which indicates continuous run. In this mode, the centrifuge will run until manually stopped. **To start a run, press the control knob.**

When the pre-selected time expires, the centrifuge will stop automatically. **To stop the centrifuge prior to the expiration of set time, press the control knob.**

The centrifuge may be operated for a short run by pressing and holding the control knob. The centrifuge will continue to run as long as the control knob is depressed and the time, in seconds will count up on the time display.

6 SERVICE AND MAINTENANCE

6.1 Centrifuge service

The brushless motor in the **PerfectSpin 24 Plus** requires no routine maintenance. Any required service should be performed by authorized, qualified personnel only. Repairs performed by unauthorized personnel may void the warranty.

6.2 Cleaning the centrifuge

Thorough cleaning not only has its purpose in hygiene but also in avoiding corrosion based on pollution. Always keep the centrifuge housing, rotor chamber, rotor and rotor accessories clean. All parts should be wiped down periodically with a soft cloth. For more thorough cleaning, use a neutral cleaning agent (pH between 6 and 8) applied with a soft cloth to avoid damaging anodized parts such as rotors. Alkaline cleaning agents (pH-value > 8) must not be used. Excessive amounts of liquid should be avoided. Liquid should not come into contact with the motor. After cleaning, ensure that all parts are dried thoroughly by hand or in a warm air cabinet (maximum temperature 50 °C).

6.3 Cleaning the rotor

The rotor should be cleaned after each use. When spinning samples containing phenol or phenol chloroform, the rotor should be cleaned immediately after use.

6.4 Disinfection of the rotor

In case of infectious material spilling into the centrifuge, the rotor and rotor chamber have to be disinfected right after the run. Rotors may be autoclaved at a maximum temperature of 121 °C. The rotor should be cleaned with a universal, neutral disinfection agent, e.g. on formalin base. A disinfection spray is most suitable in order to easily reach all difficult to access spots.

ATTENTION: Before applying any other cleaning resp. decontamination method than recommended by the manufacturer, contact the manufacturer to ensure that it will not damage the unit or the rotor.

6.5 Replacing fuses

Check the fuse when it is recommended in the Troubleshooting Guide located in this manual. The fuse holder is located in the power inlet on the rear of the unit. Disconnect the power cord from the power inlet. Open the fuse holder drawer by inserting a small screwdriver under the tab and prying it open. Remove the innermost (operative) fuse from its retaining tabs and replace the fuse if necessary. A spare fuse is located in the outermost chamber of the fuse drawer. Replace only with a fuse of exactly the same value as the original. (Fuse type may be found in the Technical data section of this manual.)

7 TROUBLESHOOTING

Please refer to this guide before calling for service.

Centrifuge will not start

Possible reason:	No power supply
Solution:	Check that power is being supplied to the outlet Check that the power cord is plugged into both the wall outlet and the back of the centrifuge Check that power cord is not damaged
Possible reason:	Blown fuse
Solution:	Check fuse and replace if necessary

Lid lock will not release

Possible reason:	Defective lid lock
Solution:	Open manually and have unit serviced
Possible reason:	No power from PC board
Solution:	Call for service
Possible reason:	Lid lock is jammed
Solution:	Call for service
Possible reason:	Centrifuge is not receiving power
Solution:	See "Centrifuge will not start"

Centrifuge cannot be started, although power is on

Possible reason:	Lid not closed correctly
Solution:	Close lid correctly
Possible reason:	No speed or time has been selected
Solution:	Set speed and/or time

bAL: Error Message Indicates imbalance

Possible reason:	Tubes not inserted symmetrically in rotor holes
Solution:	Load tubes symmetrically (see 4.4 Loading the rotor)
Possible reason:	Sample liquid in tubes not balanced
Solution:	Make sure the same volume of liquid is in each tube
Possible reason:	Defective or improperly adjusted balance sensor
Solution:	Call service

Lid: (Error message) Lid not closed

Possible reason:	Lid not closed completely
Solution:	Close lid
Possible reason:	Lid lock or sensor defective
Solution:	Call service

Other Error Messages

Er on Display

Solution:	Press time or speed knob to clear error Call service
-----------	---

8 TECHNICAL SUPPORT

Should you have any questions about the **PerfectSpin 24 Plus** or its accessories, please visit www.peqlab.com to find the respective contact person.

Please have the unit's serial number (located on the back panel of the instrument) available when calling. Should an item require return to PEQLAB for service, a decontamination form must be completed first by the user. Items sent without decontamination form will not be accepted.

All returns must be pre-approved by PEQLAB!

9 DETERMINATION OF G-VALUES

The centrifuging radius of the 1.5 / 2.0 ml rotor is 8.4 cm. The **PerfectSpin 24 Plus** has an automatic g-force conversion program, so g-values are automatically calculated and can be displayed on the control panel, for this centrifuging radius. If adapters or smaller tubes are used, the centrifuging radius changes as does the g-force. **Figure 3** can be used to manually determine g-values for any centrifuging radius.


Figure 3: Nomograph – Relative Centrifugal Force

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D PEQLAB Biotechnologie GmbH, 91052 Erlangen, Freecall (D): 0800 100 20 16, info@peqlab.de, www.peqlab.de
AT PEQLAB Biotechnologie GmbH, 6404 Polling, Tel: +43 (0) 5238 84 169, info@peqlab.at, www.peqlab.at
UK PEQLAB Ltd., Southampton SO31 7ZN, Freephone (UK): 0808 202 1302, info@peqlab.co.uk, www.peqlab.co.uk
USA PEQLAB LLC, Wilmington, DE 19810, Toll-Free (US): 877 737 5220, info@peqlab.us, www.peqlab.us